


WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

SERIES SPECIFIC

27 821 979












Product version up to 3/12/2018




- LAMINAR FLOW, max. flow rate 15 l/min (at 3 bar)
- with volume control
- rosette 60 x 60 mm
- holder for Kneipp hose
- 1/2" wall elbow
- without hose

Intrinsic protection against back flow.

Product change as of 1.4.2023: max. flow 15 l/min

	Chrome	27 821 979-00
	Brushed Platinum	27 821 979-06
	Platinum	27 821 979-08
	Dark Chrome	27 821 979-19
	Brushed Light Gold	27 821 979-27
	Brushed Durabrass (23kt Gold)	27 821 979-28
	Matte Black	27 821 979-33
	Brushed Champagne (22kt Gold)	27 821 979-46
	Champagne (22kt Gold)	27 821 979-47
	Brushed Chrome	27 821 979-93
	Brushed Dark Platinum	27 821 979-99

Recommended miscellaneous




Concealed wall elbow -

- max. recess depth 163 mm
- min. recess depth 90 mm

35 085 970 90

Required miscellaneous




WATER TUBE Kneipp hose white

2000 mm - Chrome

- with holding ring
- hose nozzle with sleeve M33

28 285 979-00

Required miscellaneous



WATER TUBE Kneipp hose black

2000 mm - Chrome


- with holding ring
- hose nozzle with sleeve M33

28 286 979-00


WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

SERIES SPECIFIC

27 821 979 Product version up to 3/12/2018
mm [inches]



Flow rate chart



Codes & Standards

EN 1717



WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

SERIES SPECIFIC

27 821 979 Product version up to 3/12/2018


Certificates and sustainability

Belgaqua_22_280_2

WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

SERIES SPECIFIC

27 821 979 Product version up to 3/12/2018



Spare parts list

No.	Item Number	Name	Quantity used	Delivery time
	8 09 24 03 151 90	nipple	1.00	2
	6 90 31 11 112 00 90	pin	1.00	2
	12 90 14 10 049 00 90	seal	1.00	2
	1 09 20 97 021-00	handle	1.00	10
	13 05 30 31 025 00 90	fixing set	1.00	2
	5 90 23 01 045 00 90	back-flow preventer	1.00	2
	14 09 27 97 022-00	rosette	1.00	10
	16 09 31 11 148 90	pin	2.00	2
	4 09 11 03 135-00	connection	1.00	10
	7 09 14 10 081 90	seal	1.00	2
	11 09 14 10 119 90	seal	1.00	2
	3 90 23 01 146 00 90	flow reducer	1.00	2
	9 09 31 11 040 90	pin	1.00	2
	15 08 17 20 726 90	holder	2.00	2
	17 09 17 20 023-00	holder	1.00	10
	10 09 27 97 021-00	rosette	1.00	10
	2 90 90 03 145 00 90	top	1.00	2