

WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

IMO

27 821 979



Fits: IMO, MEM, CL.1, CYO

- LAMINAR FLOW, max. flow rate 15 l/min (at 3 bar)
- with volume control
- rosette 60 x 60 mm
- holder for Kneipp hose
- 1/2" wall elbow
- without hose

Intrinsic protection against back flow.

Product change as of 1.4.2023: max. flow 15 l/min

| | | |
|--|-------------------------------|---------------|
| | Chrome | 27 821 979-00 |
| | Brushed Platinum | 27 821 979-06 |
| | Platinum | 27 821 979-08 |
| | Dark Chrome | 27 821 979-19 |
| | Brushed Light Gold | 27 821 979-27 |
| | Brushed Durabass (23kt Gold) | 27 821 979-28 |
| | Matte Black | 27 821 979-33 |
| | Brushed Champagne (22kt Gold) | 27 821 979-46 |
| | Champagne (22kt Gold) | 27 821 979-47 |
| | Brushed Chrome | 27 821 979-93 |
| | Brushed Dark Platinum | 27 821 979-99 |

Recommended miscellaneous



Concealed wall elbow -

- max. recess depth 163 mm
- min. recess depth 90 mm

35 085 970 90

Required miscellaneous



WATER TUBE Kneipp hose white 2000 mm - Chrome

- with holding ring
- hose nozzle with sleeve M33

28 285 979-00

Required miscellaneous



WATER TUBE Kneipp hose black 2000 mm - Chrome


- with holding ring
- hose nozzle with sleeve M33

28 286 979-00


WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

IMO

27 821 979
mm [inches]



Flow rate chart



Codes & Standards

EN 1717



WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

IMO

27 821 979

Certificates

Belgaqua_22_280_2

WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome


IMO

27 821 979

Product version up to 3/12/2018

Product version from 3/12/2018 to 3/6/2023

Product version from 3/6/2023





WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Chrome

IMO

27 821 979

Spare parts list

Product version up to 3/12/2018

Product version from 3/12/2018 to 3/6/2023

Product version from 3/6/2023

| No. | Item Number | Name | Quantity used | Delivery time |
|-----|---------------|------------|---------------|---------------|
| | 28 486 979-00 | connection | 8 | 10 |