













WATER TUBE Kneipp wall elbow with hose holder with individual rosettes -
Matte Black

SERIES SPECIFIC

27 821 979	Product version from 3/6/2023	<ul style="list-style-type: none"> • LAMINAR FLOW, max. flow rate 15 l/min (at 3 bar) • with volume control • rosette 60 x 60 mm • holder for Kneipp hose • 1/2" wall elbow • without hose <p>Intrinsic protection against back flow.</p> <p>Product change as of 1.4.2023: max. flow 15 l/min</p>	<div>  Matte Black </div> <div>  Chrome </div> <div>  Brushed Platinum </div> <div>  Platinum </div> <div>  Dark Chrome </div> <div>  Brushed Light Gold </div> <div>  Brushed Durabrass (23kt Gold) </div> <div>  Brushed Dark Brass </div> <div>  Brushed Champagne (22kt Gold) </div> <div>  Champagne (22kt Gold) </div> <div>  Brushed Chrome </div> <div>  Brushed Dark Platinum </div>	<div>27 821 979-33</div> <div>27 821 979-00</div> <div>27 821 979-06</div> <div>27 821 979-08</div> <div>27 821 979-19</div> <div>27 821 979-27</div> <div>27 821 979-28</div> <div>27 821 979-39</div> <div>27 821 979-46</div> <div>27 821 979-47</div> <div>27 821 979-93</div> <div>27 821 979-99</div>
------------	-------------------------------	--	--	---

WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Matte Black

SERIES SPECIFIC

27 821 979 Product version from 3/6/2023
mm [inches]



Flow rate chart



Codes & Standards

EN 1717



WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Matte Black

SERIES SPECIFIC

27 821 979 Product version from 3/6/2023

Certificates and sustainability


Belgaqua_22_280_2

WATER TUBE Kneipp wall elbow with hose holder with individual rosettes - Matte Black

SERIES SPECIFIC

27 821 979 Product version from 3/6/2023

Parts for other
finishes can be found
here: [Chrome](#)



Spare parts list

No.	Item Number	Name	Quantity used	Delivery time
2	90 90 03 145 00 90	top	1.00	2
15	08 17 20 726 90	holder	2.00	2
16	09 31 11 148 90	pin	2.00	2
5	90 23 01 045 00 90	back-flow preventer	1.00	2
17	09 17 20 023-33	holder	1.00	30
8	09 24 03 151 90	nipple	1.00	2
3	90 23 01 146 00 90	flow reducer	1.00	2
7	09 14 10 081 90	seal	1.00	2
4	09 11 03 135-33	connection	1.00	60
12	90 14 10 049 00 90	seal	1.00	2
10	09 27 97 021-33	rosette	1.00	60
1	09 20 97 021-33	handle	1.00	60
13	05 30 31 025 00 90	fixing set	1.00	2
6	90 31 11 112 00 90	pin	1.00	2
9	09 31 11 040 90	pin	1.00	2
11	09 14 10 119 90	seal	1.00	2
14	09 27 97 022-33	rosette	1.00	60